

# Biocard™ Quant

## Rapid Test Reader



Troponin I  
Myoglobin  
CK-MB  
CRP  
hsCRP  
Myeloperoxidase  
(MPO)

- **Operating System: Windows CE**
- **Power supply: 100V-240V AC**
- **Keyboard: Touch Screen**

Omega Diagnostics (Irl) Ltd.  
Unit 14, Bullford Business  
Campus,  
Kilcoole  
Co. Wicklow  
Tel: + 353 1 287 8158  
Fax: + 353 1 287 8149  
Email: [sales@omegairl.com](mailto:sales@omegairl.com)

# Biocard<sup>TM</sup> Quant

## Rapid Test Reader



- Quantitative results for Troponin I, Myoglobin, CK-MB, hsCRP, CRP, and MPO tests
- Desktop stand-alone instrument with embedded thermal printer
- Connectivity via USB, serial and ethernet ports
  - Patient records up to 1000 patients
  - High Precision and Accuracy
  - CE marked
- Patented dynamic laser scanning technology

Omega Diagnostics (Irl) Ltd  
Unit 14, Bullford Business  
Campus Kilcoole  
Co. Wicklow  
Ireland

Phone: +353 1 287 8158  
Fax: +353 1 287 8149  
Email: [sales@omegairl.com](mailto:sales@omegairl.com)

